

## AD-A259 050 MENTATION PAGE

Form Approved  
OMB No. 0704-0188

2



estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Avenue, Washington, DC 20540.

1. REPORT DATE 18 Dec 92		3. REPORT TYPE AND DATES COVERED Final: 30 Sept 91-30 Sept 92	
4. TITLE AND SUBTITLE Rapid detection of infectious agents with a biosensor-based nucleic acid hybridization assay		5. FUNDING NUMBERS Contract #N00014-91-C-0279	
6. AUTHOR(S) Karen M. Hahnenberger, Ph.D.		8. PERFORMING ORGANIZATION REPORT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Molecular Devices Corporation 4700 Bohannon Drive Menlo Park, CA 94025		10. SPONSORING/MONITORING AGENCY REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES) Naval Medical Research Institute Building 17A-Room 316 Bethesda, Maryland 20814-5055 Attn: Dr. Mitchell Carl		11. SUPPLEMENTARY NOTES The view, opinions and/or findings contained in this report are those of the author(s) and should not be construed as an official Department of the Army position, policy, or decision, unless so designated by other documentation.	
12a. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution unlimited.		12b. DISTRIBUTION CODE	
13. ABSTRACT (Maximum 200 words) <p>The goal of this project has been to optimize the Threshold DNA probe assay as a method for quantitating specific sequences after amplification by the polymerase chain reaction. Increasing the probe concentration 5-fold in the standard scheme B assay results in more rapid annealing of the probes to the target DNA, decreasing the overall assay time to less than 1 hour. A second assay format was developed in which the capture reagent and enzyme conjugate are added simultaneously. This scheme A format gives comparable results to those obtained using scheme B, and shortens the assay time to less than 30 minutes. Studies were also initiated to increase the detection sensitivity of the Threshold DNA probe assay by introducing multiple labels into the probe sequence by PCR. Data obtained after the reporting period officially ended indicate that 10 genomic targets can be detected after 25 cycles of PCR by using longer probes containing multiple fluorescein labels.</p> <p>92 12 23 094</p>			
14. SUBJECT TERMS DNA probes, PCR, hybridization, detection sensitivity		15. NUMBER OF PAGES 26	
17. SECURITY CLASSIFICATION OF REPORT UNCLASSIFIED		16. PRICE CODE	
18. SECURITY CLASSIFICATION OF THIS PAGE UNCLASSIFIED		19. SECURITY CLASSIFICATION OF ABSTRACT UNCLASSIFIED	
20. LIMITATION OF ABSTRACT UL			

NSN 7540-01-280-5500

Standard Form 298 (Rev. 2-89)  
Prescribed by ANSI Std. Z39-18  
298-102289  
92-32737

## GENERAL INSTRUCTIONS\*FOR COMPLETING SF 298

The Report Documentation Page (RDP) is used in announcing and cataloging reports. It is important that this information be consistent with the rest of the report, particularly the cover and title page. Instructions for filling in each block of the form follow. It is important to *stay within the lines* to meet optical scanning requirements.

**Block 1. Agency Use Only (Leave blank).**

**Block 2. Report Date.** Full publication date including day, month, and year, if available (e.g. 1 Jan 88). Must cite at least the year.

**Block 3. Type of Report and Dates Covered.** State whether report is interim, final, etc. If applicable, enter inclusive report dates (e.g. 10 Jun 87 - 30 Jun 88).

**Block 4. Title and Subtitle.** A title is taken from the part of the report that provides the most meaningful and complete information. When a report is prepared in more than one volume, repeat the primary title, add volume number, and include subtitle for the specific volume. On classified documents enter the title classification in parentheses.

**Block 5. Funding Numbers.** To include contract and grant numbers; may include program element number(s), project number(s), task number(s), and work unit number(s). Use the following labels:

C - Contract	PR - Project
G - Grant	TA - Task
PE - Program Element	WU - Work Unit Accession No.

**Block 6. Author(s).** Name(s) of person(s) responsible for writing the report, performing the research, or credited with the content of the report. If editor or compiler, this should follow the name(s).

**Block 7. Performing Organization Name(s) and Address(es).** Self-explanatory.

**Block 8. Performing Organization Report Number.** Enter the unique alphanumeric report number(s) assigned by the organization performing the report.

**Block 9. Sponsoring/Monitoring Agency Name(s) and Address(es).** Self-explanatory.

**Block 10. Sponsoring/Monitoring Agency Report Number.** (If known)

**Block 11. Supplementary Notes.** Enter information not included elsewhere such as: Prepared in cooperation with...; Trans. of...; To be published in.... When a report is revised, include a statement whether the new report supersedes or supplements the older report.

**Block 12a. Distribution/Availability Statement.** Denotes public availability or limitations. Cite any availability to the public. Enter additional limitations or special markings in all capitals (e.g. NOFORN, REL, ITAR).

DOD - See DoDD 5230.24, "Distribution Statements on Technical Documents."

DOE - See authorities.

NASA - See Handbook NHB 2200.2.

NTIS - Leave blank.

**Block 12b. Distribution Code.**

DOD - Leave blank.

DOE - Enter DOE distribution categories from the Standard Distribution for Unclassified Scientific and Technical Reports.

NASA - Leave blank.

NTIS - Leave blank.

**Block 13. Abstract.** Include a brief (*Maximum 200 words*) factual summary of the most significant information contained in the report.

**Block 14. Subject Terms.** Keywords or phrases identifying major subjects in the report.

**Block 15. Number of Pages.** Enter the total number of pages.

**Block 16. Price Code.** Enter appropriate price code (*NTIS only*).

**Blocks 17. - 19. Security Classifications.** Self-explanatory. Enter U.S. Security Classification in accordance with U.S. Security Regulations (i.e., UNCLASSIFIED). If form contains classified information, stamp classification on the top and bottom of the page.

**Block 20. Limitation of Abstract.** This block must be completed to assign a limitation to the abstract. Enter either UL (unlimited) or SAR (same as report). An entry in this block is necessary if the abstract is to be limited. If blank, the abstract is assumed to be unlimited.

**RAPID DETECTION OF INFECTIOUS AGENTS  
WITH A BIOSENSOR-BASED NUCLEIC ACID  
HYBRIDIZATION ASSAY**

Final Report  
Contract # N00014-91-C-0279  
Reporting Period 9/30/91-9/30/92

Karen M. Hahnenberger, Ph.D.  
Molecular Devices Corporation  
4700 Bohannon Drive  
Menlo Park, CA 94025

Accession For	
NTIS CRA&I	<input checked="" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution /	
Availability Codes	
Dist	Avail and / or Special
<b>A-1</b>	

DTIC QUALITY INSPECTED 1

## **ABSTRACT**

The goal of this project has been to optimize the Threshold DNA probe assay as a method for quantitating specific sequences after amplification by the polymerase chain reaction. Increasing the probe concentration 5-fold in the standard scheme B assay results in more rapid annealing of the probes to the target DNA, decreasing the overall assay time to less than 1 hour. A second assay format was developed in which the capture reagent and enzyme conjugate are added simultaneously. This scheme A format gives comparable results to those obtained using scheme B, and shortens the assay time to less than 30 min. Studies were also initiated to increase the detection sensitivity of the Threshold DNA probe assay by introducing multiple labels into the probe sequence by PCR. Data obtained after the reporting period officially ended indicate that 10 genomic targets can be detected after 25 cycles of PCR by using longer probes containing multiple fluorescein labels.

## **MATERIALS AND METHODS**

Biotinylated and single fluoresceinated probes were synthesized by Applied Biosystems Inc. The sequences of the pGEM-specific probes are given in reference 1. The anthrax-specific probes were obtained from NMRI. Dual fluoresceinated probes were synthesized by Clontech. The standard scheme B probe assay protocol (prior to the modifications described in this report) is as follows: The target DNA to be detected is diluted into a final volume of 25  $\mu$ l of hybridization buffer (30 mM sodium phosphate, 450 mM NaCl, 3 mM EDTA, 0.25% Triton X-100, pH 7.4). Probes are diluted in hybridization buffer so that 25  $\mu$ l contains  $1.8 \times 10^{11}$  molecules of each probe. 25  $\mu$ l of probe solution and 25  $\mu$ l of target solution are mixed in PCR tubes. In a thermal cycler, the reactions are denatured at 100°C then annealed at 53°C. The tubes are then transferred to ice. 100  $\mu$ l of capture reagent containing 2  $\mu$ g of streptavidin in assay buffer (10 mM sodium phosphate, 1 mM EDTA, 150 mM NaCl, 0.05% Tween 20, 0.1% BSA, pH 7.0) is added to each tube and the samples are transferred to the filtration unit wells of the Threshold instrument. Vacuum is applied at low setting (~100  $\mu$ l/min) and then 1 ml of wash buffer (10 mM sodium phosphate, 1 mM EDTA, 150 mM NaCl, 0.05% Tween 20, pH 6.5) is filtered through each well under high vacuum (~500  $\mu$ l/min). The vacuum is turned off and 200  $\mu$ l of  $\alpha$ -fluorescein-urease conjugate in assay buffer is added to each well and incubated for 30 min. The remaining conjugate is filtered under low vacuum, 1 ml of wash

buffer is filtered through under high vacuum, and the sticks are then stored in wash buffer until read.

## **RESULTS**

### **Optimization of scheme B**

The standard scheme B protocol involves a 15 minute denaturation of the target DNA to be detected and a 15 minute annealing step with the biotinylated and fluoresceinated probes (1). The annealing time could be reduced if the concentration of probes was increased. In order to determine if increasing the probe concentration would increase the background signal, various probe concentrations were tested in hybridization buffer (no target DNA present) in the scheme B format. Increasing the concentration of both probes up to 10-fold did not significantly increase the background signal (data not shown).

The effect of increasing the probe concentration on the kinetics of hybridization to  $2 \times 10^9$  targets was examined. Using the standard probe concentration ( $1.8 \times 10^{11}$  molecules of each probe per assay), at least 2-5 minutes of annealing time is required to produce maximum Threshold signal (figure 1). However, by increasing the probe concentration 5 or 10-fold, maximum Threshold signal is obtained when the samples are transferred directly to ice from the denaturation temperature, essentially reducing the annealing time to 0 minutes (figure 1). The standard scheme B assay now uses  $9 \times 10^{11}$  molecules of each probe per assay (5X the original concentration).

The effect of the length of time the  $\alpha$ -fluorescein-urease enzyme conjugate incubates with the captured probe-target complex on the Threshold stick on the signal was evaluated using the new standard amount of probes. It was found that to obtain maximum signal at the high end of the standard curve, the conjugate must incubate at least 20-30 minutes on the membrane (figure 2). At the low end of the standard curve, however, it appears that a 10 minute incubation is sufficient. The final optimized scheme B protocol is given in appendix I.

### **Development and optimization of scheme A**

A significant reduction in the total time required for the Threshold DNA probe assay would be achieved if the 30 minute incubation time with the enzyme conjugate could be eliminated. It was postulated that if the enzyme conjugate

were added to the probe-target complex in solution, rather than after the complex had been filtered through the Threshold membrane, the fluoresceins on the probe would be more accessible and the conjugate would bind much faster. A new assay format was devised in which the enzyme conjugate is added to the samples at the same time as the capture reagent.

However, in this format any fluorescein probe not annealed to target will still be present when the conjugate is added and will compete with the bound probe. An experiment was therefore done in which the probe concentration, the amount of enzyme conjugate added, and the type of conjugate used were varied in order to determine which combination would yield the highest signal in the scheme A format. Conjugates were made with polyclonal  $\alpha$ -fluorescein or monoclonal  $\alpha$ -fluorescein. The monoclonal  $\alpha$ -fluorescein antibody was obtained from Scripps Clinic for government investigational purposes. As shown in figure 3, the highest signal for  $2 \times 10^9$  target molecules in scheme A was achieved with 200  $\mu$ l of the monoclonal  $\alpha$ -fluorescein-urease conjugate in combination with  $9 \times 10^{10}$  molecules of each probe (0.5X).

The kinetics of probe hybridization with  $9 \times 10^{10}$  molecules of each probe in the scheme A format was examined with  $8 \times 10^9$  target molecules. It was found that a 5 minute annealing time is necessary in the scheme A format to get maximum signal (figure 4). The effect of the length of time the enzyme conjugate is incubated with the probe-target mixture prior to filtering was also examined. As shown in figure 5, the conjugate incubation time in scheme A does not have a significant effect on the signal in the probe assay at the low end of the standard curve. The optimized scheme A protocol is given in appendix II.

The optimized scheme B and scheme A protocols were run side by side and as shown in figure 6, the two formats give comparable results, although the signals in the scheme B format are slightly higher at the high end of the standard curve. A possible explanation for this observation is that, at high target numbers in the scheme A protocol, the second strand of the target may displace the bound probes during the 5 min annealing period, resulting in a lower signal. Since there is no annealing period in the scheme B protocol, the reduction in signal due to strand displacement is not observed.

### **Construction of model target plasmid pKHant**

The initial studies performed with the anthrax-specific probes used the 622 bp PCR product which had been gel purified and quantitated by Tim Reif and

Suresh Pillai at NMRI. In order to ensure a steady supply of easily quantitated target DNA for the experiments discussed in this report, a portion of the 622 bp anthrax PCR product was subcloned into the vector pGEM-3Z (available from Promega). The gel purified fragment was digested with *Rsa* I, which cuts within the leftward PCR primer BACA1F1, and *Bam* HI, which cuts at the 3' end of the fluorescein probe sequence. The vector pGEM-3Z was digested with *Bam* HI and *Hinc* II, both of which cut within the polylinker sequence. (The enzymes *Rsa* I and *Hinc* II leave compatible ends which can ligate to each other). The vector was also treated with calf intestinal alkaline phosphatase to prevent re-ligation of vector without insert. The prepared vector and insert were ligated and transformed into *E. coli* strain DH5 $\alpha$ . A restriction map of the resulting plasmid, pKHant, is shown in figure 7. A large scale preparation of this plasmid was made and the DNA was quantitated by OD<sub>260</sub>. The quantitated plasmid was digested with *Bam* HI and *Sph* I to release the insert and this digest was then used as the target in all future assays using the anthrax-specific probes.

### Dual fluoresceinated probes

In the standard Threshold DNA probe assay, one of the probes contains a single fluorescein capable of binding a single  $\alpha$ -fluorescein-urease conjugate. An obvious way to increase the sensitivity of the assay, then, would be to increase the number of fluoresceins on the probe which would result in more conjugate binding to the probe-target complex. To test this hypothesis, the 30-mer probe for the anthrax target (PACAFL2) was re-synthesized with 2 fluoresceins, one at the 5' end and one at the 3' end. Using this dual fluoresceinated probe to detect the anthrax target would be expected to yield twice the signal obtained with the single labeled probe. However, when the single and double labeled probes were used to generate standard curves in a side by side comparison, it was found that there was no enhancement of signal with the dual fluoresceinated probe (figure 8).

Experiments were performed to confirm the presence of two fluoresceins on the newly synthesized probe and to assess the availability of these fluoresceins to binding by  $\alpha$ -fluorescein-biotin or  $\alpha$ -fluorescein-urease. These studies indicated that there are two fluoresceins present on each probe and that they are both available for binding, but that there is significant steric hindrance in binding two  $\alpha$ -fluorescein-urease conjugates only 30 bases apart. (This was found to be true both when the probe was free and when annealed to target).

Since there are 10 base pairs per turn of double-stranded DNA it would be expected that the 2 fluoresceins on the probe would be located on the same side of the helix. Perhaps if the fluoresceins were located on opposite sides of the double helix when annealed to target, the steric hindrance in binding two conjugate molecules would be eliminated. To test this, a 35-mer probe with a fluorescein at the 5' end and at the 3' end was synthesized. This dual fluoresceinated probe, however, was also found not to enhance the signal in the assay (data not shown).

### **Reduction of background signal**

Several experiments were performed to try to lower the background signal in the probe assay by reducing the non-specific binding of the enzyme conjugate to the Threshold membrane. Threshold sticks were pretreated by filtering through 0.5 ml of different concentrations of various blocking agents before filtering the sample. These included non-fat dry milk, casein, and Denhardt's solution. Although pre-treatment with these agents did reduce the background signal, they also reduced the specific signal, and were not investigated further. However, simply filtering 0.5 ml of wash buffer through the sticks (vacuum on high setting) before applying the sample was found to decrease the background signal by 15-20% without affecting the specific signal (data not shown).

### **Introduction of multiple fluoresceins via PCR**

Recently, kits have become available from Amersham which permit the incorporation of fluorescein-dUTP into DNA during random prime labeling or PCR. Experiments were designed to determine if the sensitivity of the Threshold DNA probe assay could be increased by incorporating multiple fluoresceins into longer, double-stranded probes via PCR or by incorporating multiple fluoresceins directly into the anthrax target during the standard 25 cycles of PCR.

PCR reactions were set up using the 5'-fluoresceinated probe (PACAFL2) and oligo BACA6R1 as the primers with varying ratios of FI-dUTP:dTTP. As shown in figure 9A, lane 2, these primers amplify an approximately 114 bp fragment from the anthrax target DNA. Amplification is also quite efficient when a ratio of 75% FI-dUTP to 25% dTTP is used (figure 9A, lane 3). The PCR product is larger in this reaction, presumably due to the incorporation of multiple fluoresceins. However, using 100% FI-dUTP in PCR yields almost no product (figure 9A, lane 4).



The PCR reaction containing 75% FI-dUTP was then used as the fluoresceinated probe in the Threshold probe assay. Since the amount of fluoresceinated probe produced by PCR was not quantitated, an arbitrary amount of the PCR reaction was added, along with the standard amount of the 30-mer biotinylated probe, to the samples. It was found that the signals obtained with the PCR generated probes were higher if a protocol which is a combination of schemes A and B is used. The scheme A/B protocol is given in appendix III. Results of a scheme A/B assay using the 114 bp FI-PCR reaction compared to the single fluoresceinated 30-mer are shown in figure 9B. Use of the FI-PCR reaction as the fluorescein probe enhances the signal approximately 2-fold.

An alternative to generating multiply fluoresceinated probes is to incorporate FI-dUTP directly into the genomic PCR reactions. The Threshold DNA probe assay would then be run by adding just the biotinylated probe. PCR reactions with varying numbers of anthrax genomic targets were performed in the presence and absence of FI-dUTP (50%) for 25 cycles. Aliquots of the reactions were analyzed on a 1.2% agarose gel stained with ethidium bromide (figure 10A). The reactions containing  $10^4$  starting genomic targets clearly show a visible product on the gel (figure 10A, lanes 3 and 8). However, the reaction containing FI-dUTP gives a much fainter band which indicates that the PCR amplification is less efficient when 50% of the dTTP is replaced with FI-dUTP. The reaction containing  $10^3$  starting targets but no FI-dUTP shows a very faint band. All other reactions show no discernible product band.

A scheme B assay was run on the FI-PCR reactions using only the biotinylated 30-mer probe and the results are shown in figure 10B. The signals obtained for the genomic FI-PCR reactions are essentially the same as those obtained for the control reactions using both the fluoresceinated and biotinylated 30-mer probes. However, since the gel analysis demonstrated that the FI-PCR reactions amplified much less efficiently, there is less target in these reactions than in the controls. Also, the FI-PCR reaction containing  $10^3$  starting targets shows no band on a gel while the control reaction shows a faint band, yet the FI-PCR reaction gives a signal of 807.8 uv/sec in the probe assay while the control reaction gives a signal of 552.1 uv/sec.

Experiments were performed to try to improve the efficiency of PCR in the presence of FI-dUTP, as has been done previously for incorporation of digoxigenin-dUTP (2, 3). Reactions were set up with  $10^4$  anthrax genomic targets and varying ratios of dTTP to FI-dUTP. After 25 cycles of amplification,

aliquots were analyzed on an agarose gel (figure 11A). It was found that a 1:1 ratio of FI-dUTP to dTTP was optimal for PCR (lane 6), but the control reaction containing 100% dTTP still amplified more efficiently (lane 3). However, when the probe assay was run on these reactions, the reaction containing 50% FI-dUTP gave a signal 3 times higher than that of the control reaction (13139.0 uv/sec compared to 4423.1 uv/sec; figure 11B) even though the gel analysis indicates that the control reaction contains more target. It was also observed that the intensity of the product band decreased in reactions with less than 50% FI-dUTP (figure 11, lanes 6-8). A possible explanation is that the intensity of the product bands on a gel is a combination of ethidium bromide fluorescence and that of the incorporated fluorescein. Therefore, there may actually be more DNA in lane 8 but the combined fluorescence of the DNA and incorporated fluorescein in lane 6 makes it appear that there is more DNA present.

After many experiments, it was concluded that introduction of multiple fluoresceins directly into the genomic target during PCR is too inconsistent to be used as a routine method for increasing the sensitivity of the Threshold DNA probe assay. Experiments performed after the conclusion of the contract between NMRI and MDC indicate that a much greater and reproducible enhancement of sensitivity is achieved by generating probes containing multiple fluoresceins via PCR. These experiments will be described in a separate report.

## **CONCLUSIONS**

The standard Threshold DNA probe assay has been optimized and total assay time reduced to less than 1 hour. A second assay format has been developed which reduces the total assay time to less than 30 minutes with no loss in sensitivity. Initial studies indicate that assay sensitivity can be greatly increased by incorporation of multiple fluoresceins into probes generated by PCR.

## **REFERENCES**

1. Olson, J.D., Panfili, P.R., Zuk, R.F. & Sheldon, E.L. (1991). Quantitation of DNA hybridization in a silicon sensor-based system: application to PCR. *Molecular and Cellular Probes* 5, 351-358.
2. Stürzl, M., Oskoui, K.B. & Roth, W.K. (1992). 'Run-off' polymerization with digoxigenin labelled nucleotides creates highly sensitive and strand specific hybridization probes: synthesis and application. *Molecular and Cellular Probes* 6, 107-114.
3. An, S.F., Franklin, D. & Fleming, K.A. (1992). Generation of digoxigenin-labelled double-stranded and single-stranded probes using the polymerase chain reaction. *Molecular and Cellular Probes* 6, 193-200.

## **LEGENDS TO FIGURES**

**Figure 1. A.**  $1.8 \times 10^{11}$  (1X),  $9 \times 10^{11}$  (5X), or  $1.8 \times 10^{12}$  (10X) molecules of each probe for the anthrax model system were mixed with  $2 \times 10^9$  molecules of the 622 bp anthrax PCR product (gel purified) in a total volume of 50  $\mu$ l of hybridization buffer. Reactions were then denatured at 100°C for 5 min and then annealed at 53°C for various times (0-15 min) before transferring to ice. Samples were then processed as in scheme B.

**Figure 2.** Four sets of reactions were set up using the pGEM/HhaI target and probes. After the samples were filtered onto the Threshold sticks, conjugate was added to one of the four sticks. After 10 min conjugate was added to a second stick and 10 min later to a third. After a total of 30 min conjugate was added to the fourth stick and all sticks were filtered on low vacuum simultaneously and then washed and read as usual.

**Figure 3.** Reactions contained  $2 \times 10^9$  pGEM/HhaI targets and various concentrations of probes (0.1X-5X). After denaturation, annealing, and addition of capture reagent, different amounts of ILA conjugate or monoclonal conjugate #143 were added. The no target negative controls were set up only with 0.5X and 1X concentration of probes.

**Figure 4.** Reactions contained  $8 \times 10^9$  pGEM/HhaI targets and  $9 \times 10^{10}$  molecules of each probe (0.5X). Reactions were denatured at 100°C for 5 min and then annealed at 53°C for 0, 1, 2, 5, 10, or 15 min before transferring to ice. Capture reagent and either 500  $\mu$ l of ILA conjugate or 200  $\mu$ l of monoclonal conjugate were added. No target negative control reactions were annealed for 5 min.

**Figure 5.** Four sets of reactions were set up with the pGEM target and probes. After denaturation and annealing, capture reagent and monoclonal conjugate were added to one set of reactions. After 10 min capture reagent and conjugate were added to the second set of reactions and 10 min later to the third. After a total of 30 min capture reagent and conjugate were added to the final set of reactions. All sets were then processed simultaneously.

**Figure 6.** Standard curves were generated using the optimized scheme A and scheme B protocols (appendices I and II) side by side with the pGEM model system.

**Figure 7.** The plasmid pKHant contains a 585 bp portion of the 622 bp anthrax PCR product cloned into the vector pGEM-3Z. The location of the ampicillin resistance gene (Amp), origin of replication (ori), and lacZ gene (containing the polylinker cloning site) are shown. Also shown are representative restriction sites within the anthrax fragment.

**Figure 8.** The scheme B format was used to generate standard curves for the anthrax target (gel purified) with both the single and dual fluorescein probes. Note that although the specific signals are comparable for both probes, the background signal for the dual fluorescein probe is higher (~500 uv/sec).


**Figure 9. A.** PCR reactions (50  $\mu$ l) contained 200  $\mu$ M dNTPs, 0.5  $\mu$ M of each primer,  $10^7$  targets (gel purified 622 bp fragment), and 2.5 units of Pfu DNA polymerase (Stratagene; buffer supplied by manufacturer). Fluorescein-11-dUTP (1mM solution) was purchased from Amersham. Cycling parameters were as follows: 25 cycles of 94°C, 1 min; 55°C, 1 min; 72°C, 2 min and 1 cycle of 72°C, 7 min. 10  $\mu$ l of each reaction was analyzed on 1.2% Seakem GTG agarose (FMC BioProducts) in 1X Tris Acetate EDTA buffer (25X stock purchased from AMRESCO). **B.** The scheme A/B format was used to generate standard curves using the single fluorescein probe and 10  $\mu$ l of the FI-PCR reaction shown in lane 3 above. pKHant digested with *Bam* HI and *Sph* I was used as the target.

**Figure 10. A.** PCR reactions (50  $\mu$ l) contained 200  $\mu$ M dNTPs, 0.5  $\mu$ M of primers (BACA1F1 and BACA6R1), 5 units of Taq polymerase (Promega; buffer supplied by manufacturer). The target was 0, 10, 100, 1000, or 10,000 copies of the *Bacillus anthracis* genome (supplied by Suresh Pillai at NMRI). FI-PCR reactions contained 100  $\mu$ M dTTP and 100  $\mu$ M FI-dUTP. Cycling conditions were as follows: 1 cycle of 94°C 5 min; 25 cycles of 94°C 1 min, 55°C 1 min, 72°C 2 min; 1 cycle of 72°C 7 min. 10  $\mu$ l aliquots of each reaction were analyzed on 1.2% agarose in 1X TAE buffer. **B.** The scheme B format was used to assay 10  $\mu$ l aliquots of each PCR reaction. Both biotin and fluorescein probes were added to the control reactions while only the biotin probe was added to assay the


reactions containing FI-dUTP.

**Figure 11. A.** PCR reactions were set up with  $10^4$  *Bacillus anthracis* genomic targets and varying ratios of FI-dUTP to dTTP. Reactions were cycled and analyzed as described in the legend to figure 11A. **B.** 5  $\mu$ l aliquots of each reaction were assayed using the scheme B format using only the biotin probe. The positive control (dTTP only) and negative control (no target) PCR reactions were assayed with both the fluorescein and biotin probes.


**Figure 1**



**Scheme B: Effect of conjugate incubation time**







**Figure 3**

**Figure 4**





Figure 3



Optimized Scheme A vs. Scheme B





Figure 7



Dual FI probe vs. single FI probe




Figure 9

A.



Lane

Sample

1

123 bp ladder

2

Control (100% dTTP)

3


75%FI-dUTP:25%dTTP

4


100% FI-dUTP

B.

Single FI probe vs. 114bp FI-PCR probe



**A.**




**Lane**

**Sample**

1	123 bp ladder
2	No target
3	$10^4$ targets
4	$10^3$ targets
5	$10^2$ targets
6	$10^1$ targets
7	No targets + FI-dUTP
8	$10^4$ targets + FI-dUTP
9	$10^3$ targets + FI-dUTP
10	$10^2$ targets + FI-dUTP
11	$10^1$ targets + FI-dUTP


**B.**

**Anthrax genomic PCR reactions +/- FI-dUTP**






A.



Lane	Sample
1	123 bp ladder
2	No target
3	Control: 200 $\mu$ M dTTP
4	80%FI-dUTP:20%dTTP
5	67%FI-dUTP:33%dTTP
6	50%FI-dUTP:50%dTTP
7	33%FI-dUTP:67%dTTP
8	20%FI-dUTP:80%dTTP

B.

Effect of FI-dUTP:dTTP ratio on anthrax genomic PCR efficiency



## APPENDIX I

### SCHEME B DNA PROBE ASSAY

- |  |         |
|--|---------|
| 1. Denature target at 100°C  | 5 min   |
| 2. Anneal probes ( $9 \times 10^{11}$ molecules of each probe) at 53°C | 0 min*  |
| 3. Add 100µl of 10X capture reagent, filter, and wash                  | ~10 min |
| 4. Incubate with 200µl αFl-urease conjugate (ILA)                      | 30 min  |
| 5. Filter and wash   | ~5 min  |
| 6. Read sticks   | ~5 min  |

\* In the optimized scheme B protocol the reactions are transferred to ice directly from the denaturation temperature and then immediately processed.

## APPENDIX II

### OPTIMIZED SCHEME A PROBE ASSAY

- |    |   |         |
|----|---|---------|
| 1. | Denature target at 100°C  | 5 min   |
| 2. | Anneal probes ( $9 \times 10^{10}$ molecules of each) at 53°C                                     | 5 min   |
| 3. | Add 100µl of 10X capture reagent and 200µl of αFl-urease conjugate (monoclonal), filter, and wash | ~10 min |
| 4. | Read sticks   | ~5 min  |

### APPENDIX III

#### SCHEME A/B DNA PROBE ASSAY

- |  |         |
|--|---------|
| 1. Denature target at 100°C  | 5 min   |
| 2. Anneal probes ( $9 \times 10^{10}$ molecules of biotin probe and 10 $\mu$ l of FI-PCR reaction) at 53°C | 5 min   |
| 3. Add 100 $\mu$ l of 10X capture reagent, filter, and wash  | ~10 min |
| 4. Incubate with 200 $\mu$ l $\alpha$ FI-urease conjugate (ILA)  | 30 min  |
| 5. Filter and wash   | ~5 min  |
| 6. Read sticks   | ~5 min  |